



AF-6 Medium

Freshwater medium adapted for Botryococcus

Reference: Provasoli, L. and Pinter, I. J. (1959). Artificial media for fresh-water algae: problems and suggestions. In: Tryon, C. A., Jr. and Hartmann, R. T. (eds.) The Ecology of Algae. Spec. Pub. No. 2, Pymatuning Laboratory of Field Biology, Pittsburgh pp 84-96.

STOCK SOLUTIONS	CONCENTRATION: g L ⁻¹ DEIONISED WATER (dH ₂ O)	QUANTITY FOR MEDIUM
1. MES	<i>add reagent directly to medium</i>	400 mg
2. NaNO ₃	140 g	1.0 mL
3. NH ₄ NO ₃	20 g	1.1 mL
4. MgSO ₄ ·7H ₂ O	80 g	0.375 mL
5. KH ₂ PO ₄	10 g	1.0 mL
6. K ₂ HPO ₄	5 g	1.0 mL
7. CaCl ₂ ·2H ₂ O	29.4 g	0.34 mL
8. Fe citrate	<i>see recipe below</i>	0.222 mL
9. Biotin	0.1 g	20 µL
10. Thiamine HCl	0.1 g	0.1 mL
11. Vitamin B ₆	0.1 g	10 µL
12. Vitamin B ₁₂	0.1 g	10 µL
13. P IV Metal Mix	<i>see recipe below</i>	5.0 mL

Store all stock solutions in the refrigerator.

Fe citrate solution

Add both constituents to 1 L of dH₂O and autoclave to dissolve. Store solution in the dark.

CONSTITUENT	QUANTITY
Ferric citrate	9.0 g
Citric acid	9.0 g

P IV Metal Mix

Add EDTA to 1 L of dH₂O and dissolve, then add other constituents to dissolve and autoclave to sterilise. Store solution in the dark.

CONSTITUENT	QUANTITY
Na ₂ EDTA.2H ₂ O	1.0 g
FeCl ₃ .6H ₂ O	196 mg
MnCl ₂ .4H ₂ O	36.0 mg
ZnSO ₄ .7H ₂ O	22.0 mg
CoCl ₂ .6H ₂ O	4.0 mg
Na ₂ MoO ₄ .2H ₂ O	2.5 mg

To prepare AF-6 Medium (1 L)

- Add MES (1) to 990.0 mL dH₂O
- Add each stock solution (2 – 13) in the stated volume to the MES solution.
- Adjust the pH to 6.5 – 7 with 20% NaOH.
- Autoclave at 121°C (15 psi for 20 mins).

CONTACT US

t 1300 363 400
+61 3 9545 2176
e csiroenquiries@csiro.au
w www.csiro.au

For further information

Australian National Algae Culture Collection
w www.csiro.au/en/Research/Collections/ANACC

Ian Jameson
Director
t +61 3 6232 5117
e ian.jameson@csiro.au

**Australian National Algae
Supply Service**

Cathy Johnston
Manager
t +61 3 6232 5316
e cathy.johnston@csiro.au