

## Embedding for light microscopy of GUS-stained tissues

Use this protocol for best structural preservation

1. Immediately after X-Glc stain, proceed as for standard LRW embedding, starting with 3-5 rinses in 25 mM phosphate buffer, pH 7.2, 15 min each rinse.
2. Fix in 3% glutaraldehyde in 25 mM phosphate buffer, pH 7.2, at room temperature for 90 min or longer. If necessary (should not be necessary after X-Glc stain), vacuum infiltrate first for a few min until material sinks.
3. Rinse in 25 mM phosphate buffer 4 times over approx 60 min.
4. Dehydrate in ethanol series:

Concentration of ethanol	Dehydration times
25%	1-2 h
50%	1-2 h
70%	overnight at 4°C
95%	2 h
100%	1-2 h, repeat at least 3 times

5. Infiltrate in LR White resin (medium grade):

Concentration of resin	Infiltration times
2:1 ethanol : LRW	8 h to overnight
1:2 ethanol : LRW	8 h to overnight
100% LRW	At least overnight, then repeat 3 times, at least 8 h each

5. Embed in rubber moulds or gelatine capsules. In each case, a thin layer can be polymerized first to facilitate positioning of tissue in mould. Alternatively, for large tissues, flat embed in foil planchettes. Polymerise a thin layer of resin on the bottom of the planchettes before embedding (at 70°C under nitrogen gas).
6. Polymerise in vacuum oven at 70°C for about 90 minutes under nitrogen gas (polymerisation is complete when resin is hard), adjust flow rate to give 250-300 kPa on the nitrogen tank gauge.

Waste LRW and LRW-solvent mixtures go into glass jar in fume hood.

Waste LRW-contaminated pipettes go into plastic container in fume hood.

### **Embedding for light microscopy of GUS-stained tissues**

Use this protocol if tissue has already been rinsed and cleared in 70% ethanol

1. Proceed as for standard LRW embedding:

Dehydrate in ethanol series:

<b>Concentration of ethanol</b>	<b>Dehydration times</b>
70%	overnight at 4°C
95%	2 h
100%	1-2 h, repeat at least 3 times

2. Infiltrate in LR White resin (medium grade):

<b>Concentration of resin</b>	<b>Infiltration times</b>
2:1 ethanol : LRW	8 h to overnight
1:2 ethanol : LRW	8 h to overnight
100% LRW	At least overnight, then repeat 3 times, at least 8 h each

5. Embed in rubber moulds or gelatine capsules. In each case, a thin layer can be polymerized first to facilitate positioning of tissue in mould. Alternatively, for large tissues, flat embed in foil planchettes. Polymerise a thin layer of resin on the bottom of the planchettes before embedding (at 70°C under nitrogen gas).
6. Polymerise in vacuum oven at 70°C for about 90 minutes under nitrogen gas (polymerisation is complete when resin is hard), adjust flow rate to give 250-300 kPa on the nitrogen tank gauge.

Waste LRW and LRW-solvent mixtures go into glass jar in fume hood.

Waste LRW-contaminated pipettes go into plastic container in fume hood.

**Notes:**

- i. Larger and denser tissues need longer in every stage. For example, cereal grains can be fixed overnight or longer, need long vacuum infiltration of fixative, need at least 8 hours in each dehydration step (preferably longer), and need slow infiltration with resin – smaller increases in resin concentration and longer at each stage – at least a day (24 h) in each. Very hard/dense tissues may need several days to several weeks for resin infiltration.
- ii. Larger tissues (as above) may need to be vacuum infiltrated in resin to remove any remaining solvent.
- iii. Impermeable tissues need slower and longer infiltration or the tissues will shrink and be deformed – impermeable walls allow solvent (EtOH) to diffuse out much more rapidly than the resin monomers can diffuse in, hence cell collapse. For example, *Arabidopsis* roots show cell collapse if the resin concentration is increased too rapidly – commonly seen in GUS-stained tissues if not careful.
- iv. Hard tissues such as dry grains may also need to be embedded in harder resin – hard grade rather than medium grade.