

Guide to the See3 Viewing/Scoring/Optimisation tool in C3

Janet Newman

March 2017

Introduction:

See3 is the software used by the Collaborative Crystallisation Centre (C3) in Melbourne to view and to score images. The software also allows users to create reports and design optimisations. The software is based on the “XtalTrak” application developed by Rigaku, but we have worked extensively on it to add functionality (and get rid of some niggles). Ideas for further improvements always welcome!

Geek stuff:

The program is written in the ASP MVC .NET Framework using KendoUI for the user interface. Alex Kruger, Nick Rosa and Marko Ristic have all worked to transform the original Rigaku code into See3.

About the Login page:


Our Twitter feed is on the login page, this is where you will find the latest info from C3 – look here if something doesn’t seem to be working: if we haven’t tweeted about it, we probably don’t know about it – please let us know (email c3@csiro.au) about any problems you are having.

Using the software:


Navigate to See3.csiro.au in the browser of your choice, and log on with your C3 username and password.



This will open up the “Plates” tab, which is dominated by a grid on the left hand side of the page. The grid can take several seconds to load. The bar to the right of the grid may be moved by dragging the middle of the bar (circled in red, below) with the left mouse button depressed. Many of the panes can be resized the same way.



The grid can be modified as you wish – the changes will only persist for the current session. You can change the order of the grid columns, sort by a column, group by a column, or filter a column for a particular string (filtering shown below).



Clicking on a row in the grid will open up the latest inspection for that plate. There are a number of options on this view:

There are three tabs on the upper, left hand side (circled in red, below) – ‘Plates’, ‘Optimisation’ and ‘Clipboard’. Call these *VIEWS*

There are three tabs in the upper, middle screen (circled in green) ‘General’, ‘Plate’, ‘Drop’. Call these *TABS*

There are four pulldown menus on the upper, right hand side (circled in purple) ‘Sort by’, ‘Score’, ‘Type’ ‘Subwells’

There are 6 tabs in the middle of the screen, below the drops (circled in navy blue) ‘Inspections’, ‘Conditions’, ‘Scoring’, ‘History’, ‘Tools’ and ‘Reports’. Call these *OPTIONS*

See3 v4.0.3.20192

Catalase: MC103211 (13: AL1 / 20' / Score / PlateName: 35474_shotgun_20_1

Plates Optimization Clipboard

General Plate Drop

Sort By: Position Current Score Mono 2mm All

Barco...	T (°C)	User	Project	Sample	Screen	Last Inspected	Date
MC103214	20	sumeet.bal...	Cytochrome C	Ferritin	shotgun	03/11/2017 04:22	01/13
MC103213	20	sumeet.bal...	Cytochrome C	Cytochrome...	shotgun	03/11/2017 04:12	01/13
MC103212	20	sumeet.bal...	Catalase	Irinin	shotgun	03/11/2017 03:57	01/13
MC103211	20	sumeet.bal...	Catalase	Catalase	shotgun	03/11/2017 03:47	01/13
MC103210	20	sumeet.bal...	Ribonuclease A	Thermolysin	shotgun	03/11/2017 03:37	01/13
MC103209	20	sumeet.bal...	Ribonuclease A	Ribonuclease...	shotgun	03/11/2017 03:26	01/13
MC103204	20	sumeet.bal...	Thaumatin	Thaumatin	shotgun	03/11/2017 03:15	01/06
MC103203	20	sumeet.bal...	concanavalin A	alpha-Lactal...	shotgun	03/11/2017 01:43	01/06
MC103202	20	sumeet.bal...	concanavalin A	Concanaval...	shotgun	03/11/2017 01:33	01/06
G0000181	20	sumeet.bal...	Projects	AtzB	OPFF_shotgun_sumee...	03/02/2017 01:07	12/20
MC103166	20	sumeet.bal...	AtzB	AtzB	OPFF_shotgun_sumee...	03/01/2017 00:14	12/20
G0000180	20	sumeet.bal...	Projects	Cysteine Hy...	OPFF_shotgun_sumee...	02/15/2017 08:27	12/07
MC103101	20	sumeet.bal...	cysteine hydrolase	Cysteine Hy...	OPFF_shotgun_sumee...	02/15/2017 07:14	12/07
MC103067	20	sumeet.bal...	h-Carbonic Anhydra...	h-Carbonic ...	OPFF_shotgun_sumee...	02/03/2017 19:30	11/10
G0000179	20	sumeet.bal...	Projects	h-Carbonic ...	OPFF_shotgun_sumee...	02/03/2017 19:28	11/10
G0000177	20	sumeet.bal...	Projects	Myoglobin	OPFF_shotgun_sumee...	01/26/2017 17:39	11/10
MC102990	20	sumeet.bal...	myoglobin	Myoglobin	OPFF_shotgun_sumee...	01/26/2017 15:29	11/10
MC102949	20	sumeet.bal...	h-Carbonic Anhydra...	h-Carbonic ...	shotgun	01/19/2017 11:20	11/10
MC102947	20	sumeet.bal...	AtzB	Cysteine Hy...	shotgun	01/18/2017 19:02	11/09
MC102946	20	sumeet.bal...	AtzB	AtzB	shotgun	01/18/2017 15:32	11/09
MC102941	20	sumeet.bal...	myoglobin	Myoglobin	shotgun	01/16/2017 15:11	11/07
G0000174	20	sumeet.bal...	Projects	Ins	OPFF_shotgun_sumee...	01/13/2017 20:03	11/04

Inspections Conditions Scoring History Tools Reports

Inspection	Type	Date	Elapsed	S	Expects	Views	Scores
13	Auto	03/11/2017	36d	1	1	0	0
12	Auto	03/03/2017	49d	1	0	0	0
11	Auto	02/24/2017	41d	1	0	0	0
10	Auto	02/17/2017	34d	2	0	0	0
9	Auto	02/10/2017	27d	1	2	0	0
8	Auto	02/03/2017	21d	1	0	0	0
7	Auto	01/27/2017	14d	1	4	0	0

Plates VIEW, Plate TAB:

Sort by

'Sort by' allows you to arrange the images either by position in the plate, or by an assigned score. By default, human scores are shown first, and better scores are shown first. Below is the same plate as above, but with the 'Sort by Score' option selected:

Catalase: MC103211 (13: AL1 / 20' / Score / PlateName: 35474_shotgun_20_1

Welcome Janet Newman Logout

General Plate Drop


Sort By: Score Current Score Mono 2mm All

Score


In the C3 database, inspections rather than drops are scored. This gives us the ability in See3 to choose between displaying the 'Current Score' (i.e., the most recent score associated with the inspection being viewed) and the 'Last Score', which is the most recent score associated with any inspection of that drop. By default, the 'Last Score' is shown. If there are both machine and human scores for a droplet, one can view only human scores by selecting the 'Last User Score' option.

Type

In C3 we generally collect zoomed colour images collected with visible light. Our imaging schedules include other types of images – for example, images collected with UV light, or with polarising filters in place. If there is more than one type of image collected at the same time (i.e. in the same inspection), they can be accessed by choosing the appropriate image type from the 'Type' pulldown. The example below shows that there were both visible and UV images for a particular inspection, and the visible option has been selected.




The first inspection of a plate will show both an overview inspection and the zoomed inspection.




Subwells

C3 uses two subwell plates (SwissSci SD-2 plates, aka 'MRC plates'); by default all drops are shown. If you want to only view the drops in the top (or the bottom) subwell then choose the appropriate subwell set from this pulldown



Plates VIEW, Drop TAB:

In this arrangement, the 'Sort' menu changes to 'Split', which allows the user to see two images from the same inspection (for example, Visible and UV images) side by side.



Most of the six *OPTIONS* give information or functions the Drop *TAB*. In the Drop *TAB*, the image can be zoomed in using the mouse wheel, and the green scale bar will adjust according to the zoom

level. Once zoomed, holding down the left mouse button allows you to drag the drop around. The zooming and re-centring on the image persist until another plate is loaded from the plate grid.

Plates *VIEW*, Drop *TAB*,

Inspections OPTION

Shows the list of inspections for a selected plate, and gives the Inspection Number, the Date of Inspection, the time since the plate was set up, the number of sub-inspections, how many times an inspection has been viewed, and if any scores were associated with an inspection. This is the default *OPTION*.

Inspection	Type	Date	Elapsed	Subinspects	Views	Scores
14	Auto	03/11/2017	63d	1	0	0
13	Auto	03/03/2017	56d	1	0	0
12	Auto	02/24/2017	48d	1	0	0
11	Auto	02/17/2017	41d	1	1	0
10	Auto	02/10/2017	34d	2	5	1
9	Auto	02/03/2017	28d	1	0	0
8	Auto	01/27/2017	21d	1	8	6

Conditions OPTION

Shows details about the chemistry of the droplet (what is in the reservoir), information about the sample (concentration, sample_ID), the ratio of the sample to the reservoir and information about additives, if appropriate. By default, the treeview for the reservoir contents is open, to see more information about the sample, click on the small arrow to the left of the sample (circled in red, below)

Inspection	Conditions	Scoring	History	Tools	Reports
Recipe	Type	Source	Volume (µL)		
	Sample	Concanavalin A	0.20		
	Crystallant	shotgun	0.10		
	Conc	Units	Name *	pH	Group
	20	w/v	poly(ethylene glycol) 3000		Polymer
	0.1	M	trisodium citrate-citric acid	5.5	Buffer

Click the Recipe button (circled green, above) to find out how to make up 1 mL of reservoir solution using C3 stocks. Use this recipe when booking a 'One Tube' request on the TECAN

Recipe
>
21.0 ul of citric acid "pH 2.3" (1M)
400.0 ul of poly(ethylene glycol) 3000 (50% w/v)
79.0 ul of (tri)sodium citrate "pH 8.1" (1M)
500.0 ul of H2O


Scoring OPTION

There are 16 possible user scores one can associate with a droplet – see the table below. The order of the scores in the table below shows the ordering of scores when the 'Sort by Score' menu option is chosen. The descriptions associated with a score are somewhat arbitrary, as one project's 'Crystals***' might be another project's 'Crystal*', but the idea is that the more stars there are, the better. When a score is associated with a droplet, a coloured rim is drawn around the crystal image. The colours are grouped – yellow/red colours are associated with crystals, blue colours with precipitate, green colours for other results.

SCORE	DESCRIPTION	SHORTCUT
SHOOT ME	In plate screening marker	m


CRYSTALS ***	Worth harvesting	6
CRYSTALS **	Needs some work, but very promising	5
CRYSTALS *	too many, too small, needs a lot of optimising	4
CRYSTALLINE	Microcrystalline, optimise only if desperate	c
SPHERULITES	Things that look like meteors	8
OTHER	Unknown	9
SALT XTAL?	Putative salt crystals	s
GOOD PRECIPITATE	Grainy, shiny, birefringent	1
PRECIPITATE	All other precipitates	3
BAD PRECIPITATE	Brown, smoky, denatured	2
PHASE SEPARATION	Oily, immiscible	7
CLEAR	Clear droplet	0
CLEAR WITH STUFF	Clear drops with dust, fibres, or other non-crystalline matter	j
CONFIRMED SALT	Confirmed by UV or diffraction	b
NULL EXPERIMENT	This experiment didn't work (dried, not put down well)	n

One can assign a score to a droplet by clicking on one of the coloured boxes in the *Scoring OPTIONS* (circled red below), or by using the keyboard shortcut given in the table. When a score is associated with a drop, the viewing software will advance to the next image, and the scorer and score date will show up in the Scoring *OPTION*.




History *OPTION*

Shows the time-course of the droplet, that is, all the images associated with the droplet to date. Clicking on one of the thumbnails in the History *OPTION* will move it to the main viewing pane.



Tools OPTION


The measure tool is selected by clicking on the ruler icon (circled red, below left), and is de-selected by clicking on the four-arrow icon (circled green, below right). This option allows one to measure an object in a drop. The green ruler lines on the image persist until another plate is loaded from the plate grid.



The re-centre/re-size tool is selected by clicking on the crosshair icon (circled in purple, below)




Once clicked, use the mouse to re-centre/re-zoom the droplet (some guesswork might be involved, as is the case above and below).



Click the green tick or the red cross to exit the re-centre tool. If you selected the tick the new centre and extent will be used for the next inspection of that drop, and the icon will appear in the drop image to confirm that the drop has been re-centred.


The forward, back and play buttons (circled below) move you through the images – the delay between the images can be set here as well.



The autoplay function can also be toggled on and off with the space bar, and the keyboard arrow keys will also move you between images.

Reports *OPTION*

There are two reports, one which gives an overview of the selected plates, and one which gives a more detailed report of images that have been placed on the clipboard. Click on the “Generate” button opposite the report that you desire.




The Overview Report gives a .pdf report of the plate as it is laid out



The Hit Report gives details about individual droplets that have been selected. At the end of the Hit Report, a summary of the chemistry found in the drops selected for the report is given, an example is shown below.

See3 Hit Report

MC103207 A1.1 Inspection #14



Sample Name: Lyszyme
20 mg/ml S. lyszyme

Screen Name: Lyszyme QA Screen
DIO: n/a C sodium HEPES (pH 7)
O.S. New/v C polyethylene glycol 8000
1.1 M C sodium citrate-citric acid (pH 7)


Volumes (ul)
Sample: 015
Crystallant: 015

Project (path): Lyszyme\CLG\Alicia.a1
Project (date, time): 15/01/2017 15:03/2017
Elapsed Time: 3.3 Days

Scores: (type, user, date, inspection #)

Visible image: FOV 2.5mm

MC103207 A2.1 Inspection #14



Sample Name: Lyszyme
20 mg/ml S. lyszyme

Screen Name: Lyszyme QA Screen
DIO: n/a C lithium chloride
20 New/v C polyethylene glycol 8000
1.0 M C sodium citrate-citric acid (pH 4)


Volumes (ul)
Sample: 015
Crystallant: 015

Project (path): Lyszyme\CLG\Alicia.a1
Project (date, time): 15/01/2017 15:03/2017
Elapsed Time: 3.3 Days

Scores: (type, user, date, inspection #)

Visible image: FOV 2.5mm

MC103207 A3.1 Inspection #14



Sample Name: Lyszyme
20 mg/ml S. lyszyme

Screen Name: Lyszyme QA Screen
DIO: n/a C magnesium chloride
0.1 M C tris(hydroxymethyl)aminomethane (pH 8.5)
20 New/v C polyethylene glycol 8000

Volumes (ul)
Sample: 015
Crystallant: 015

Project (path): Lyszyme\CLG\Alicia.a1
Project (date, time): 15/01/2017 15:03/2017
Elapsed Time: 3.3 Days


Scores: (type, user, date, inspection #)

Visible image: FOV 2.5mm

CSRD 1 Generated by Janet Newman on March 15, 2017

See3 Hit Report

MC103207 A4.1 Inspection #14



Sample Name: Lyszyme
20 mg/ml S. lyszyme

Screen Name: Lyszyme QA Screen
DIO: n/a C sodium chloride
0.1 M C sodium acetate-acetic acid (pH 4.5)

Volumes (ul)
Sample: 015
Crystallant: 015

Project (path): Lyszyme\CLG\Alicia.a1
Project (date, time): 15/01/2017 15:03/2017
Elapsed Time: 3.3 Days

Scores: (type, user, date, inspection #)

Visible image: FOV 2.5mm

Chemistry Range Table

Class	Avg Conc	Conc Range	Units	Name	Count
Buffer	100		mM	sodium HEPES (pH 7)	1
Buffer	0.1		M	tris-chloride (pH 8.5)	1
Buffer	0.1		M	sodium acetate-acetic acid (pH 4.5)	1
Other	0.5		w/v	polyethylene glycol 8000	1
Salt	1.1		M	sodium malonate-malonic acid (pH 7)	1
Salt	1		M	lithium chloride	1
Salt	0.2		M	magnesium chloride	1
Precipitant	20		w/v	polyethylene glycol 8000	1
Precipitant	20		w/v	polyethylene glycol 8000	1
Precipitant	1		M	sodium chloride	1
Fundamental Mixtures	100		mM	trisodium citrate-citric acid (pH 4)	1

CSRD 2 Generated by Janet Newman on March 15, 2017

Clipboard VIEW

The clipboard allows the user to assemble a set of images (droplets) that are interesting – these can be from one or more plates. The clipboard is cleared each time a user logs out of the system. A drop is added to the clipboard by holding down the CTRL key whilst clicking the appropriate image. This works both from the Plate or the Drop VIEW.

See3 v1.0.3.30332 Lyszyme: MC103207 (14), A1.1 / 20 / Score / PlateName: 35471_Lyszyme QA Screen_20_1 Welcome Janet Newman [Logout](#)


Plates Optimization Clipboard

Barcode	Well	Inspection
MC103207	A1.1	14

(ctrl+click images to add them to the clipboard)

[Clear Clipboard](#) [Compare first two drops](#)

General Plate Drop [Split](#) Current Score Mono 2mm All



Inspections Conditions Scoring History Tools Reports

Inspection	Type	Date	Elapsed	Subimages	Views	Scores
14	Auto	03/15/2017	62d	1	0	0
13	Auto	03/08/2017	55d	1	0	0
12	Auto	03/01/2017	48d	1	0	0
11	Auto	02/22/2017	41d	1	1	0
10	Auto	02/15/2017	34d	1	0	0

Once the appropriate set of drops has been assembled, these can be used in either optimisation or to generate a Hit report. One can add all the images from the History *OPTION* to the clipboard, then use the Hit Report to get a record of the time course of a droplet

The screenshot displays the See3 software interface. On the left is a data table with columns: Barcode, Well, Inspection, and a Remove button. The table lists 8 rows of data for MC103207 in well A1.1, with inspection counts ranging from 7 to 14. Below the table are buttons for 'Clear Clipboard' and 'Compare first two drops'. On the right is a large circular microscope image showing a petri dish with several dark spots. Below the image is a timeline of smaller images with dates from 10/03/2017 to 7/21/2017. A red box highlights the 'Add all to clipboard' button in the bottom left of the image area.

Barcode	Well	Inspection	Remove
MC103207	A1.1	14	Remove
MC103207	A1.1	13	Remove
MC103207	A1.1	12	Remove
MC103207	A1.1	11	Remove
MC103207	A1.1	10	Remove
MC103207	A1.1	9	Remove
MC103207	A1.1	8	Remove
MC103207	A1.1	7	Remove

(ctrl+click images to add them to the clipboard)

Clear Clipboard Compare first two drops

General Plate Drop Split Current Score Mono 20mm All

Scoring History Tools Reports

Add all to clipboard

10/03/2017 12: 1/03/2017 11: 22/03/2017 10: 15/03/2017 9: 8/02/2017 8: 1/03/2017 7: 25/01/2017 6: 21/01/2017

Optimisation VIEW

The Optimisation VIEW has three areas, the **Information Pane**, found on the left hand side of the VIEW, the **Factor Grid**, on the upper right hand side of the VIEW and the **Plate Grid**, on the lower right side of the VIEW.

The screenshot shows the See3 v1.0.3.30332 interface. The **Information Pane** on the left includes tabs for Plates, Optimization, and Clipboard. Under Optimization, there are buttons for Randomize, Add Factor Group, and Add Chemical. Below these are fields for Plate Size (96), Additive (Select), Dilution (0), and a checkbox for Include Initial Conditions. The Factor Grid on the right lists three factor groups: Primary (polyethylene glycol 8000), Salt (magnesium chloride), and Buffer (tris chloride). Each group has a table of parameters including Name, Avg Conc, Start Conc, End Conc, Units, Avg pH, Start pH, End pH, Weight, Vary, Steps, and Distribution. The Plate Grid at the bottom is a 12x8 grid of wells, each containing a horizontal bar chart with red and green segments, representing the current instantiation of the Factor Grid.

The **Information Pane** reads from and writes to the underlying database. This pane allows one to select the 'Plate Size' - the size of the optimisation: 24 well (4 rows x 6 columns); 48 well (6 rows x 8 columns) or most commonly 96 well (8 rows x 12 columns). Following standard convention, rows are lettered and columns are numbered and is always shown in portrait orientation with well A1 in the upper left hand corner. Thus a 96-well plate has rows A-H, columns 1-12.

In the **Information Pane** one can populate the **Factor Grid** by adding Factor Groups and Chemicals to the optimisation - manually or from the clipboard, select designs to be used as additive screens, set the Design name and save the Design and Factor Grid.

The **Factor Grid** encapsulates the concepts of the optimisation, whereas the **Plate Grid** shows the current instantiation generated by that particular **Factor Grid**. Thus the same **Factor Grid** will generate different **Plate Grids** if the 'Randomize' button on the **Information Pane** is clicked. If the 'Save Design' button on the **Information Pane** is pressed, the conditions currently shown in the **Plate Grid** will be saved. A previously saved **Factor Grid** can be reloaded by clicking the 'Load Factors' button – this will display a list of previously saved grids, which can be selected.

The **Plate Grid** gives a visual representation of the current instance of the **Factor Grid**, where the colour of the bar indicates from which Factor Group it comes, and the length of the bar gives an indication of the concentration or the pH (the distance from the lower limit). Mousing over a well will display the contents of the well. By default, any conditions used to define the optimisation will be transferred unchanged to the optimisation, and will be found in the wells starting at position C3. These wells have a pale blue background, and the starting condition is included in the information provided in the mouse over. Uncheck the 'Include Initial Conditions' checkbox in the **Information Pane** to turn off this feature (Red Arrow, below)

Randomize
Add Factor Group
Add Chemical

Plate Size:

Additive:

Dilution (%):

Include Initial Conditions: ☒

User Name:

Design Name:

Save Design
Save Factors
Load Factors
Generate from Clipboard

Name	Avg Conc	Start Conc	End Conc	Units	Avg pH
polyethylene glycol 8000					
20	16	22	w/v		

Salt

Name	Avg Conc	Start Conc	End Conc	Units	Avg pH
magnesium chloride					
0.2	0.02	0.3	M		

Buffer

Name	Avg Conc	Start Conc	End Conc	Units	Avg pH
tris chloride					
0.1	0.1	0.1	M		8.5

	1	2	3	C3 (MC103239 A3)	6
A					
B					
C					

Mousing over a particular coloured bar will highlight the other wells which contain that same chemical (black dashed border), and will highlight where that chemical is found in the Factor Grid. Conversely, mousing over a chemical in the **Factor Grid** will highlight where that chemical is found in the **Plate Grid**.

Welcome Janet Newman
Logout

Factor Group	Color	Vary By	Location	Steps	
Primary		Gaussian Random	Random	96 / 96	Delete
Name	Avg Conc	Start Conc	End Conc	Units	Avg pH
polyethylene glycol 8000					
20	16	22	w/v		
Salt		Gaussian Random	Random	96 / 96	Delete
Name	Avg Conc	Start Conc	End Conc	Units	Avg pH
magnesium chloride					
0.2	0.02	0.3	M		
magnesium acetate	0.2	0.02	0.3	M	
Buffer		Gaussian Random	Random	96 / 96	Delete
Name	Avg Conc	Start Conc	End Conc	Units	Avg pH
tris chloride					
0.1	0.1	0.1	M		8.5

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												

The **Factor Grid** consists of one or more Factor Groups. Each Factor Group contains one or more Chemical Factors. A Chemical Factor is a chemical with an associated concentration, unit and (potentially) pH value. Thus 'ammonium sulfate' is a chemical, and '1 M ammonium sulfate' is a Chemical Factor. The **Factor Grid** has a Header Bar (marked with a red arrow, below) which contains columns relevant to all Factor Groups. Each individual Factor Group has a secondary Header (green arrows) which contains columns relevant to the Chemical Factors within the Factor Group. The coloured 'Delete' buttons remove the entire Factor Group, whereas the blue 'Delete' buttons remove individual Chemical Factors from within a Factor Group.

Welcome Janet Newman

Logout

Factor Group					Color	Vary By				Location	Steps	
Primary						Gaussian Random				Random	96 / 96	<div>Delete</div>
Name	Avg Conc	Start Conc	End Conc	Units	Avg pH	Start pH	End pH	Weight	Vary	Steps	Distribution	
polyethylene glycol 8000	20	16	22	w/v	7	7	7	0.5	Conc	96	48	<div>Delete</div>
sodium malonate-malonic acid	1.1	0.88	1.21	M				0.5	Conc	96	48	<div>Delete</div>
Salt						Gaussian Random				Random	48 / 96	<div>Delete</div>
Name	Avg Conc	Start Conc	End Conc	Units	Avg pH	Start pH	End pH	Weight	Vary	Steps	Distribution	
magnesium chloride	0.2	0.02	0.3	M				0.5	Conc	96	48	<div>Delete</div>
Buffer						Gaussian Random				Random	96 / 96	<div>Delete</div>
Name	Avg Conc	Start Conc	End Conc	Units	Avg pH	Start pH	End pH	Weight	Vary	Steps	Distribution	
tris chloride	0.1	0.1	0.1	M	8.5	7.1	9.1	0.5	pH	96	48	<div>Delete</div>
sodium HEPES	100	100	100	mM	7	6.55	8.55	0.5	pH	96	48	<div>Delete</div>
Other						Gaussian Random				Random	48 / 96	<div>Delete</div>
Name	Avg Conc	Start Conc	End Conc	Units	Avg pH	Start pH	End pH	Weight	Vary	Steps	Distribution	
jeffamine ED-2001	0.5	0.05	0.75	w/v	7	7	7	0.5	Conc	96	48	<div>Delete</div>


</

Within a Factor Group one can vary the colour, the way in which concentration varies ('concentration' will be used throughout, the same logic would apply to pH values if the pH rather than concentration is being varied), and the way in which the values are arrayed in the **Plate Grid**. All factors in the group will use the same 'Colour', 'Vary By' and 'Location' options.


The 'Vary By' pulldown menu shows the different options for changing the concentration of the factors within the group. By default, the concentration is varied randomly between the upper (starting x 110%) and lower (starting x 80%) limits, but with a bias towards the starting value. This is called 'Gaussian Random'. An unbiased random selection of concentration between the upper and lower limits can be chosen with the 'Vary By' option 'Uniform Random' – these are shown as the red and blue images below.




Both the Gaussian and Uniform random selections can be ordered, so that the smallest concentration value is found in the upper right hand corner of the selection, and the largest concentration is found in the lower right. Notice that the values don't change, but the arrangement does.




'Vary By' column sets up a grid from left to right across the columns of the plate (green, below) and 'Vary By' row sets up a grid from top to bottom across the rows of the plate (yellow, below)




'Vary By' block gives four values spread over the quadrants of the plate (cyan, below), and 'Vary By' Page gives two values spread over the two halves of the plate (magenta, below)




'Vary By' fixed keeps the value of the factor constant (brown, below)




When there is more than one Chemical Factor in the Factor Group then the 'Location' option becomes relevant. This allows one to position the Chemical Factors on the plate, irrespective of how the concentration is being varied. By default, the distribution of the Chemical Factors on the plate is 'Random' (red, below – look to see which wells have the black dotted border). By selecting 'Column' under 'Location', the plate is sorted so that all the wells containing a particular Chemical Factor are co-located (blue, below – look to see which wells have the black dotted border). The distribution of the Chemical Factor across the plate is determined by the weight associated with the Chemical Factor. In the following example, each of the two Chemical Factors within the Factor group have an equal weighting of 0.5, so so each gets half the number of available columns. The order of the Chemical Factors within the group dictates where on the plate the clustering occurs – the first Chemical Factor in the Factor Group will be to the left (Location = Column or Plate) to the top (Location = Row), to the upper left (Location = Block). One can drag and drop Chemical Factors within a Factor Group to fine tune the location on the plate.




If the number of Chemical Factors or weighting of the Chemical Factors in a Factor Group does not divide evenly into the location (for example, below – three Chemical Factors of equal weight cannot be spread exactly into a Location = Row), the software will round as best it can to whole columns, rows, blocks or pages – and this will be shown in the 'Distribution' column of the Secondary Header – circled in red, below.



If starting from scratch use the 'Add Chemical' button to populate the Factor Grid – click on the 'Add Chemical' button, and a dialog will pop open




Starting from the top, populate the fields – the 'Name' field shows a list chemicals available in C3. The pulldown menu list can be filtered by including a string in the 'Name' Field. Notice that the full chemical name is used, not common abbreviations. The string 'poly' typed into the 'Name' field will return all the polyethylene glycols, but the string 'PEG' will not.




Once the chemical has been selected, chose the desired factor group (or just use the default 'New Group' – the name can be changed later) and the desired concentration. By default the limits of the chemical will be 80-100% of the given concentration UNLESS the chemical concentration is used at 0.2 M or less, in which case the default limits are set to 10-150% of the given concentration. If the pH is varied, then by default the limits are set to the [nearest pKa] +/- 1 pH unit. If these limits are not what you want, then the limits can be changed

Use the 'New Factor Group' button to create one or more Factor Groups in the Factor Grid, rename as required, then these new Factor Groups can be selected as the destination for the Chemical Factor.

A Design is saved into the database by pressing the 'Save Design' button on the **Information Pane**. When 'Save Design' is clicked, a check is made to ensure that all the stocks are available in C3 to make the required design. If the design has been made where there is a mismatch between the units of the required Chemical Factor and the units of the stock available in C3 this will cause an error. This happens when, say, 20 w/v polyethylene glycol 200 has been used as a Chemical Factor, and the stock in C3 is 80 v/v polyethylene glycol 200. More commonly, there may be 'overflow' errors, where the stocks available in C3 are not concentrated enough to be able to make the combination of chemical factors. For example, if a well was designed to contain 30 w/v polyethylene glycol 4000, 2M ammonium sulfate – this can't be made with C3 stocks of 50 w/v polyethylene glycol 4000 and 3.5M ammonium sulfate (as $30/50 + 3/3.5 > 1$). If there are less than 20 overflow errors, the software will give the option of autofixing (Click 'Autofix'), which it does by reducing the concentration of the primary factor until there is no error, if possible. Clicking 'OK' will return you to the Optimization VIEW to fix this manually.



The software will warn, but not error, if incompatible chemicals are co-located in wells. Combinations that are incompatible include magnesium salts with phosphate salts and zinc salts with sodium HEPES buffer. Press 'OK' to go back into the Optimization VIEW, or 'Ignore' to continue writing the design into the database.



The **Factor Grid** can be saved by clicking the 'Save Factors' button on the **Information Pane**. The saved **Factor Grid** will have the 'Design Name' as shown in the **Information Pane** at the time the button is clicked.